



 **MEDIREX**

Positive Patient ID
& Process Solutions

Laminated ID Bands for
Laser Printers

The Laminated Laser Wristband

BARCODE READY

The Laminated Laser Wristbands are ideal for environments that require a laminated seal to preserve patient information against fading and smudging. The self-laminated design with tamper-evident adhesive closure allows for easy application. These wristbands are compatible with laser and ink jet printers, and suitable for all hospital visits up-to 7 days. The laminated laser wristbands are available in adult and infant wristband sizes, a variety of laser label quantities on the same sheet, and various colour selections.

Patient Safety

- Laminated seal to resist moisture, alcohol, and hand sanitizer
- Hypoallergenic, latex-free, phthalate-free, and single-use disposable

Durability

- Developed with a tamper-evident adhesive closure
- Recommended for visits lasting up-to 7 days

Printer Compatibility


- Suited for all laser printer manufacturers
- Compatible with ink jet printers

Easy Application

- Wristband application instructions highlighting a 4-step self-laminating process for easy application

Highly Versatile

- Available in adult and infant wristband sizes, a variety of laser label quantities on the same sheet, and colours suited for your unique needs



Laminated Laser Adult Wristband with 8-Labels

BARCODE READY

The Laminated Laser Adult Wristband with 8-Labels is the ideal solution for low acuity visits or short stays by reducing the wastage of labels when a high label count is not required. The laminated seal preserves patient information against fading and smudging, while the tamper-evident adhesive closure is tested and proven to be durable for hospital visits up-to 7 days. A 4-step self-laminating process provides easy wristband application and trouble-free operator experience. These wristbands are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

Item No.	Colour
7505	<input type="radio"/> Solid White
7505R	<input checked="" type="radio"/> Red Border

Dimensions (L x H)

Overall Size: 279MM X 147MM [11" X 5.8"]
 Wristband: 267MM X 20MM [10.50" X 0.75"]
 Label (8): 64MM X 25MM [2.5" X 1.0"]

Packaging String: 1CASE/4PACKAGES/250EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

FURTADO, SARA
 DOB: 06/26/86 AGE: 35Y
 Unit #: 12345678

MEDIREX
 The Laminated Wristband - 7505
www.medirex.com | orders@medirex.com

STEP 1 Print ID band and labels on laser printer.
STEP 2 (See below)
STEP 3 Turn sheet over and lift ID band off sheet; start from the right and pull to the left.
STEP 4 Fold form over imaged area to create ID band.

SELF LAMINATING WRISTBAND

7505

Push lower edge of ID band down to begin separating imaged band from sheet.

Licensed under one or more of the following U.S. Patent 6,653,472, and LaserBand U.S. Patents 6,510,634,881; 6,067,739; 6,000,160; and others pending in the U.S. and Canada

Laminated Laser Family Wristbands with 1-ID Card

BARCODE READY

The Laminated Laser Family Wristbands with 1-ID Card is recommended for Labour & Delivery. This family set of wristbands vastly improves patient safety by having both adult and infant wristbands on the same sheet for printing. The laminated seal preserves patient information against fading and smudging, while the tamper-evident adhesive closure is tested and proven to be durable for hospital visits up-to 7 days. A 4-step self-laminating process provides easy wristband application and trouble-free operator experience. These wristbands are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

Item No. 7501P
Colour Solid White

Dimensions (L x H)

Overall Size: 279MM X 216MM [11.0" X 8.5"]
 Adult Wristband (2): 267MM X 20MM [10.50" X 0.75"]
 Infant Wristband (2): 137MM X 13MM [5.625" X 0.50"]
 ID Card (1): 82MM X 51MM [3.25" X 2.0"]

Packaging String: 1CASE/4PACKAGES/250EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).



Custom Products by request

At Medirex, we always get excited about the next challenge. We look forward to putting our years of experience, industry knowledge, and creativity to good use. Medirex is proud to develop personalized patient identification solutions that suit the unique needs of your hospital environment, are best-in-class, value-driven, and focused on safety.

For personalized solutions of wristband colours or label sizes and quantities, please reach out to our sales team:

✉ sales@medirex.com

☎ 416.363.9313
1.800.387.9848

For general inquiries please contact:

✉ info@medirex.com

Other ways to get in touch:

🌐 www.Medirex.com

📘 @medirexsystems

🌐 <https://ca.linkedin.com/company/medirex>

🐦 @medirexsys

